

# X3D WebGL based Volume Rendering

Luis Kabongo [lkabongo@vicomtech.org](mailto:lkabongo@vicomtech.org)


John Congote [jcongote@vicomtech.org](mailto:jcongote@vicomtech.org)

Esther Novo [enovo@vicomtech.org](mailto:enovo@vicomtech.org)



# Volume Visualization 2011

**Real-Time Interactive Visualization of Volumetric Data with WebGL**



Rotation with the mouse, and Alt+Mouse to Zoom

<http://demos.vicomtech.org/volren/index.html>



# 2012

- Implementation to support advanced medical visualization on the Web without plugins using existing standards ...



# MEDX3D

Extension of the X3D ISO standard to support advanced medical visualization functionality

```
<?xml version="1.0" encoding="utf-8"?>
<X3D version="3.2" profile="MedicalInterchange">
<Scene>
<Group>
<NavigationInfo type='EXAMINE' WALK' />
<Viewpoint description="first" position="0 0 0.6"/>
<Background skyColor="1 1 1"/>
<Transform scale=".001 .001 .001">
<!-- The volume -->
<VolumeData DEF="volume" rayStep="0.01" useSlicing="false"
  dimensions="250 250 156">
<Image3DTexture containerField="voxels" url='data/S76280/
  IM00010">
</VolumeData>
</Transform>
</Group>
</Scene>
</X3D>
```



# X3DOM (Simple Example)

Framework for integrating and manipulating X3D scenes as HTML5/DOM elements

```
<!DOCTYPE html PUBLIC "-//W3C//DTD XHTML 1.0 Strict//EN" "http://www.w3.org/TR/xhtml1/DTD/xhtml1-strict.dtd">
<html xmlns="http://www.w3.org/1999/xhtml">
    <head>
        <meta http-equiv="X-UA-Compatible" content="chrome=1" />
        <meta http-equiv="Content-Type" content="text/html; charset=utf-8" />
        <title>Hello World</title>
        <link rel="stylesheet" type="text/css" href="x3dom.css" />
        <script type="text/javascript" src="x3dom.js"></script>
    </head>
<body>
    <h1>XHTML Hello World</h1>
    <p>With X3D-namespace and case sensitive element/node names. Works with self-closing tags (e.g. Viewpoint and Material)</p>
    <X3D xmlns="http://www.web3d.org/specifications/x3d-namespace" showStat="false" showLog="false" x="0px" y="0px" width="400px" height="400px" altImg="helloX3D-alt.png">
        <Scene>
            <Viewpoint position='0 0 10' />
            <Shape>
                <Appearance>
                    <Material diffuseColor='0.603 0.894 0.909' />
                </Appearance>
                <Box DEF='box' />
            </Shape>
        </Scene>
    </X3D>
</body>
</html>
```



# MEDX3DOM(Work in progress)

```
<!DOCTYPE html PUBLIC "-//W3C//DTD XHTML 1.0 Strict//EN" "http://www.w3.org/TR/xhtml1/DTD/xhtml1-strict.dtd">

<html xmlns="http://www.w3.org/1999/xhtml">
<head>
<meta http-equiv="X-UA-Compatible" content="chrome=1" />
<meta http-equiv="Content-Type" content="text/html; charset=utf-8" />
<title>Hello World</title>
<link rel="stylesheet" type="text/css" href="x3dom.css" />
<script type="text/javascript" src="x3dom.js"></script>
</head>
<body>
    <X3D xmlns='http://www.web3d.org/specifications/x3d-namespace' showStat='true' showLog='true' width='500px' height='500px'>
        <Scene>
            <Background skyColor='0.0 0.3 0.65' />
            <Viewpoint description='Default' zNear='0.0001' zFar='100' />
            <Transform>
                <VolumeData id='volume' dimensions='4.0 4.0 4.0'>
                    <ImageTextureAtlas containerField='voxels' url='media/volume/aorta4096.png' number_of_slices='96' slices_over_x='10' slices_over_y='10' />
                    <OpacityMapVolumeStyle lightFactor='0.01'>
                        <ImageTexture containerField='transferFunction' url='transfer.png' />
                    </OpacityMapVolumeStyle>
                </VolumeData>
            </Transform>
        </Scene>
    </X3D>
</body>
</html>
```



# Our contribution in MedX3DOM

## Generation of nodes for two components

- Texturing3D

- X3DTexture3DNode
- ComposedTexture3D
- ImageTexture3D
- PixelTexture3D
- TextureCoordinate3D
- TextureCoordinate4D
- TextureTransformMatrix3D
- TextureTransform3D
- ImageTextureAtlas

Implemented

Defined

Not implemented

Abstract

- VolumeRendering


- X3DComposableVolumeRenderStyleNode
- X3DVolumeDataNode
- X3DVolumeRenderStyleNode
- BlendedVolumeStyle
- BoundaryEnhancementVolumeStyle
- CartoonVolumeStyle
- ComposedVolumeStyle
- EdgeEnhancementVolumeStyle
- IsoSurfaceVolumeData
- OpacityMapVolumeStyle
- ProjectionVolumeStyle
- SegmentedVolumeData
- ShadedVolumeStyle
- SilhouetteEnhancementVolumeStyle
- ToneMappedVolumeStyle
- VolumeData



web|3D  
CONSORTIUM

# ImageTextureAtlas for an Atlas Image Type

- Atlas
  - Composed image
  - Parameters
    - Number of slices
    - Number of rows
    - Number of columns
- Brain Atlas Image



*The texture atlas is the mosaic of all the slices of a volume in one image, the order of the images is given by the rows and columns.*



# New Node ImageTextureAtlas

```
ImageTextureAtlas: ImageTexture {  
    SFInt32 [in]      numberOfSlices      0  
    SFInt32 [in]      slicesOverX        0  
    SFInt32 [in]      slicesOverY        0  
}
```

ImageTextureAtlas node defines a texture map by specifying a single image file that contains slides of complete 3D data. Parameters indicate the information necessary for mapping voxels to the geometry.



web|3D

CONSORTIUM

# Modified OpacityMapVolumeStyle

```
OpacityMapVolumeStyle : X3DComposableVolumeRenderStyleNode {  
    SFBool [in,out] enabled          TRUE  
    SFNode [in,out] metadata         NULL  [X3DMetadataObject]  
    SFNode [in,out] transferFunction NULL  [X3DTexture2DNode,X3DTexture3DNode]  
    SFFloat [in,out] opacityFactor   0.01  [0,1]  
    SFFloat [in,out] lightFactor     0.3   [0,1]  
}
```

**OpacityFactor:** Value of contrast in each voxel of the volume


**LightFactor:** Value of radiance of light in each voxel of the volume



web|3D  
CONSORTIUM

# Brain Example

```
<X3D xmlns='http://www.web3d.org/specifications/x3d-
  namespace'
  showStat='true' showLog='true' width='500px'
  height='500px'>
<Scene>
  <Background skyColor='0.0 0.0 0.0' />
  <Viewpoint description='Default' zNear='0.0001'
  zFar='100' />
  <Transform>
    <VolumeData id='volume' dimensions='4.0 4.0 4.0'>
      <ImageTextureAtlas containerField='voxels'
      url='media/volume/brain-at_4096.jpg'
      number_of_slices='226' slices_over_x='16'
      slices_over_y='16' />
      <OpacityMapVolumeStyle>
      </OpacityMapVolumeStyle>
    </VolumeData>
  </Transform>
</Scene>
</X3D>
```




[http://demos.vicomtech.org/x3dom/test/functional/volrenOpacityTest\\_brain.xhtml](http://demos.vicomtech.org/x3dom/test/functional/volrenOpacityTest_brain.xhtml)



web|3D  
CONSORTIUM

# Skull Example

```
<X3D xmlns='http://www.web3d.org/specifications/x3d-
namespace'
showStat='true' showLog='true' width='500px'
height='500px'>
<Scene>
  <Background skyColor='0.0 0.0 0.0' />
  <Viewpoint description='Default' zNear='0.0001'
zFar='100' />
  <Transform>
    <VolumeData id='volume' dimensions='4.0 4.0
4.0'>
      <ImageTextureAtlas containerField='voxels'
url='media/volume/skull-at.png'
numberOfSlices='256' slicesOverX='16'
slicesOverY='16' />
      <OpacityMapVolumeStyle>
        </OpacityMapVolumeStyle>
    </VolumeData>
  </Transform>
</Scene>
</X3D>
```



[http://demos.vicomtech.org/x3dom/test/  
functional/volrenOpacityTest\\_skull.xhtml](http://demos.vicomtech.org/x3dom/test/functional/volrenOpacityTest_skull.xhtml)



web|3D  
CONSORTIUM

# Aorta Example

```
<X3D xmlns='http://www.web3d.org/specifications/x3d-
namespace'
showStat='true' showLog='true' width='500px'
height='500px'>
<Scene>
  <Background skyColor='0.0 0.0 0.0' />
  <Viewpoint description='Default' zNear='0.0001'
zFar='100' />
  <Transform>
    <VolumeData id='volume' dimensions='4.0 4.0
4.0'>
      <ImageTextureAtlas containerField='voxels'
url='media/volume/aorta4096.png'
numberOfSlices='96' slicesOverX='10'
slicesOverY='10' />
      <OpacityMapVolumeStyle>
      </OpacityMapVolumeStyle>
    </VolumeData>
  </Transform>
</Scene>
</X3D>
```



[http://demos.vicomtech.org/x3dom/test/functional/volrenOpacityTest\\_aorta.xhtml](http://demos.vicomtech.org/x3dom/test/functional/volrenOpacityTest_aorta.xhtml)




**web|3D**  
CONSORTIUM

# Aorta with color transfer function example

<http://demos.vicomtech.org> click on “Aorta with Transfer Function” in MEDX3DOM Section

[http://demos.vicomtech.org/x3dom/test/functional/volrenOpacityTestTF\\_aorta.xhtml](http://demos.vicomtech.org/x3dom/test/functional/volrenOpacityTestTF_aorta.xhtml)



```
INFO: Time for setup and init of GL element no. 0: 12397 ms.  
INFO: activate Bindable null  
INFO: register bindable Fog  
INFO: create new Fog for X3DFogNode-stack  
WARNING: (32,16): warning X3206: implicit truncation of vector type  
WARNING: (107,18): warning X3206: implicit truncation of vector type  
(24,5): warning X3553: can't use gradient instructions in loops with  
break, forcing loop to unroll  
INFO: activate Bindable null  
INFO: activate first X3DViewpointNode for X3DViewpointNode-stack  
INFO: activate Bindable null  
INFO: activate first X3DBackgroundNode for X3DBackgroundNode-stack  
WARNING: (34,16): warning X3206: implicit truncation of vector type  
WARNING: (36,16): warning X3206: implicit truncation of vector type  
WARNING: (34,17): warning X3206: implicit truncation of vector type
```



web|3D  
CONSORTIUM

# DICOM, NRRD, RAW Format??

- Javascript is very slow to process binary data
- The loaders of this kind of volume data are based in C and the conversion or reimplementation is not straightforward
- Other types of volume formats exists, like MHD, NetCDF, ...., impossible to implement all
- **MEDX3DOM** is for the web, using web formats





# MEDICAL IMAGING

- DICOM format is the facto standard for medical imaging, but internally is a container of several kind of data, sound, 2D images, 3D images, metadata, etc
  - [http://demos.vicomtech.org/medica/  
WtMetalImageSimpleViewer2D.wt](http://demos.vicomtech.org/medica/WtMetalImageSimpleViewer2D.wt)
  - [http://demos.vicomtech.org/medica/  
WtWebGLSimpleImageViewer.wt](http://demos.vicomtech.org/medica/WtWebGLSimpleImageViewer.wt)
- WADO: It is a webservice which provides medical images through the web retrieving from a PACS server
- MEDX3DOM: Will support WADO (Web Access to DICOM Objects)



# MEDX3DOM Future

- Work in progress
  - MPR Style
    - [http://localhost/x3domMPR/test/functional/volrenMPR\\_aorta.xhtml](http://localhost/x3domMPR/test/functional/volrenMPR_aorta.xhtml)
  - Combination with mesh models
    - [http://volumerc.org/demos/brainviewer/webgl/brain\\_viewer/brain\\_viewer.html](http://volumerc.org/demos/brainviewer/webgl/brain_viewer/brain_viewer.html)



web|3D  
CONSORTIUM


# MEDX3DOM Future

- Next steps
  - Different composition styles (MIP, X-Ray, Composed ...)
  - Lighting (Phong, Global illumination ...)
  - Animation (4D timesteps, video, flow animation ...)
    - <http://localhost/MedX3DOMGit2/test/functional/volrenTest.xhtml>
  - Data transfer optimization (streaming, compression,...)



# Beyond medical imaging

- Geovisualization
- Volume Attention Maps
- Meteorology
- Flow Visualization
- Confocal microscopy
- ....



# Keep updated!

- Our work:
  - <http://demos.vicomtech.org>
  - <http://www.volumerc.org>
- The technologies:
  - <http://www.khronos.org/webgl/>
  - <http://www.x3dom.org/>

Thank you! Questions?

